PROTHESE DISCALE / SANS CIMENT
 DISC PROSTHESIS/CEMENTLESS
 DISKUSPROTHESE / ZEMENTFREI
 PRÓTESIS DISCAL / SIN CEMENTO
 DISKUSPROTESE / CEMENTFRI
 VÁLILEVYPROTEESI / ILMAN SEMENTIÄ
 PROTEZA DISKA / BEZ CEMENTA
 PROTESI DISCALE/SENZA CEMENTO
 SKIVEPROTESE / UTEN SEMENT
 DISCUSPROTHESE / ZONDER CEMENT
 PROTEZA KRAZKA MIĘDZYKRĘGOWEGO / BEZCEMENTOWA
 PROTESE DISCAL / SEM CIMENTO
 PROTEZA DISCALĂ / FĂRĂ CIMENT
 DISKPROTES / CEMENTFRI
 PLATNICKOVA PROTÉZA/BEZCEMENTOVÁ
 DISK PROTEZA / BREZCEMENTNA
 DISK PROTEZI / CIMENTOSUZ
 THIẾT BỊ ĐĨA ĐỆM/ KHÔNG CÓ XI MĂNG
 ΠΡΟΘΕΣΗ ΔΙΣΚΟΥ / ΧΩΡΙΣ ΤΣΙΜΕΝΤΟ
 人工椎间盘假体 / 非骨水泥型
 قرص تيففري اصطناعي / من دون مادة مُلصقة

NOTICE D'INSTRUCTIONS
 INSTRUCTIONS FOR USE
 GEBRAUCHSANWEISUNG
 INSTRUCCIONES DE USO
 BRUGSANVISNING
 KÄYTTÖOHJEET
 UPUTE ZA UPORABU
 ISTRUZIONI PER L'USO
 BRUKERVEILEDNING
 GEBRUIKSAANWIJZING
 INSTRUKCJA STOSOWANIA
 FOLHETO DE INSTRUÇÕES
 MANUAL DE UTILIZARE
 BRUKSANVISNING
 NÁVOD NA POUŽITIE
 NAVODILA
 KULLANIM KILAVUZU
 HƯỚNG DẪN SỬ DỤNG
 ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
 说明书
 تعليمات الاستعمال

CP. **ESP**

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WARNINGS AND INDICATIONS FOR USE
PLEASE READ THIS DOCUMENT CAREFULLY

Do not hesitate to consult our sales department if you need further information. This document is not exhaustive, nor is it a surgical operating manual presenting details about implantation. In all cases, please consult the included documents.

INTENDED USE:

The devices covered in these instructions for use are intended for cervical spine implantation between C3 and C7 during cervical disc arthroplasty.

DESCRIPTION:

SPINE INNOVATIONS supplies a CP-ESP® cervical disc prosthesis.

The CP-ESP® cervical disc prosthesis is composed of two titanium alloy plates, separated by a deformable portion that is intended to allow the disc to move. The prosthesis has spikes on the external surface and an hydroxyapatite (HA) coating over a layer of pure titanium (T40). The deformable portion consists of an elastomeric cushion made from polycarbonate urethane.

The CP-ESP® cervical disc prosthesis is supplied sterile (sterilised using gamma irradiation) and is single-use.

The CP-ESP® cervical disc prosthesis is broadly for use at one spinal level.

References	Type	Anteroposterior length (mm)	Mediolateral length (mm)	Height (mm)
264363	CP-ESP® CEMENTLESS DISC PROSTHESIS SIZE: S1 (13x15) H5	13	15	5
264364	CP-ESP® CEMENTLESS DISC PROSTHESIS SIZE: S1 (13x15) H6	13	15	6
264365	CP-ESP® CEMENTLESS DISC PROSTHESIS SIZE: S1 (13x15) H7	13	15	7
264366	CP-ESP® CEMENTLESS DISC PROSTHESIS SIZE: S2 (14x17) H5	14	17	5
264367	CP-ESP® CEMENTLESS DISC PROSTHESIS SIZE: S2 (14x17) H6	14	17	6
264368	CP-ESP® CEMENTLESS DISC PROSTHESIS SIZE: S2 (14x17) H7	14	17	7
264369	CP-ESP® CEMENTLESS DISC PROSTHESIS SIZE: S3 (15x20) H5	15	20	5
264370	CP-ESP® CEMENTLESS DISC PROSTHESIS SIZE: S3 (15x20) H6	15	20	6
264371	CP-ESP® CEMENTLESS DISC PROSTHESIS SIZE: S3 (15x20) H7	15	20	7

IMPLANT MATERIALS:

The materials that the implants are made from are noted on each product label.

The implant is radio-opaque, except for the polycarbonate urethane.

Quantitative composition of the implanted materials in mass percentage:

References	Ti6Al4V	PCU	T40	HA
264363	81%<m<90%	5%<m<12%	3%<m<4%	2%<m<3%
264364				
264365				
264366				
264367				
264368				
264369				
264370				
264371				

Ti6Al4V: Titanium, aluminium-6 and vanadium-4 alloy

PCU: Poly Carbonate Urethane

T40: Titanium

HA: Hydroxyapatite

COMPATIBILITY BETWEEN DEVICES:

This is a monobloc device and does not require additional fixation elements.

PERFORMANCE/CLINICAL BENEFIT:

The CP-ESP prosthesis leads to reduced pain and restores the intervertebral space and the functions of the native disc. The summary of safety and clinical performance (SSCP) for the CP-ESP® cervical disc prosthesis can be found in the European database of medical devices (Eudamed):

<https://ec.europa.eu/tools/eudamed>.

This summary can be found on the Eudamed website by entering the UDI-ID for the CP-ESP® disc prosthesis: **376036923CP01A2**.

INDICATIONS FOR USE:

The CP-ESP® cervical disc prosthesis is designed for surgical treatment of degenerative disc disease (DDD) of the cervical spine in skeletally mature patients who have not responded after at least 6 months of conservative treatment. The assessment of the extent of disc involvement is based on the combination of data from the clinical examination and/or the symptoms and imaging.

The clinical examination and symptoms are defined by pain and/or a functional/neurological deficit in the neck or arm, with at least one of the following:

- Herniated disc
- Spondylarthrosis (defined by the presence of osteophytes)
- Radiculopathy
- Spinal cord compression.

The CP-ESP® cervical disc prosthesis is a medical device for primary treatment.

TARGET POPULATION:

Skeletally mature patients failing to respond to at least 6 months of conservative treatment and presenting the indications listed above and none of the contraindications listed below.

CONTRAINDICATIONS :

Specific contraindications

- Fractures, infections, tumours
- Stenosis of the spinal canal resulting from hypertrophic spondylarthrititis
- Degeneration of facet joints
- Pathological segmental instability
- Ossification of posterior longitudinal ligament

General contraindications

- Osteoporosis, osteochondrosis or severe osteopenia
- Chronic or acute local, spinal, or systemic infections
- Metabolic and systemic diseases
- Pathologies and surgical situations that preclude any benefit of spinal surgery
- Sensitization to foreign bodies resulting in reaction to implant materials
- Dependence on drugs, drug addiction or alcoholism
- Pregnancy
- Lack of cooperation from the patient

PRECAUTIONS:

A - GENERAL PRECAUTIONS

- The lifespan of implants is affected by numerous biological, biomechanical, and other factors. As a result, carefully following the indications, contraindications, precautions, and warnings for this product plays an essential role in its use. The outcome of disc prosthesis depends on the patient's medical history. Psychological preparation is indispensable. Patients must be informed of the limitations of the device, including among others, the impact of overburdening caused by weight or excessive activities. They must be advised as to how to modify their activities accordingly. Joint replacement can never reproduce the functions formerly performed by a normal, healthy joint. The patient should not hold on to unrealistic hopes about its functionality and should consult the surgeon in the event of abnormal problems experienced at the device location.
- Before clinical use, the surgeon and the surgical team must be trained in using the device and its ancillary equipment. They must have fully understood all aspects of the surgery, as well as the limitations of the device.
- Care should be taken to ensure that the components and all polished or coated surfaces are protected against abrasion, scratches, or any other damaging effects caused by abrasive or metallic objects.
- The contact zones of the implant must be clean and free of debris (bone, cement) before implantation. The surgeon should make certain to obtain even contact of the implant on the bone.
- It is advisable to make sure that the presence of other devices does not compromise the device's functional integrity.
- Implants supplied by SPINE INNOVATIONS must not be mixed with implants of another origin, or with different systems.

B - CEMENTLESS IMPLANTS

During handling, great care should be given to the coating. Any contact with the coating before implantation must be avoided. Biologically fixed implants (identification on the product label: T40 + HA) must not be cemented.

IMPLANTATION TECHNIQUE:

The specific instrument set is supplied non-sterile and is necessary to fit implants designed by SPINE INNOVATIONS. The surgical techniques are available from SPINE INNOVATIONS or its distributors.

Before the prosthesis is mounted on the introducer, adjust the selector so that the graduation marks show at least 7 mm. Open the double sterile packaging of the implant and mount the CP-ESP® disc prosthesis on the introducer. Maintaining

pressure on the lateral blades of the introducer, rotate the selector clockwise to the zero-graduation mark (the anterior edge of the prosthesis meets the endstop of the introducer).

Check the orientation of the implant before fitting it (curved surface of the disc upwards).

CONDITIONS OF USE:

The prostheses must be implanted in an operating facility, under aseptic conditions, and in compliance with hygienic practice, by a surgeon specialised in the treatment of spinal pathologies who regularly practices this type of implantation. Implantation must be carried out using the appropriate ancillary equipment, with undamaged instruments, following the recommended surgical technique and the instructions for use, to treat a patient who presents with the indications defined above.

WARNINGS:

Patients receiving an implant must be informed that the implant's lifespan may depend on their weight and activity level. Damaged or defective implants must not be used. The design of this device does not allow for its reuse, as this would risk altering its mechanical characteristics and this may impair performance. Reuse of a single-use device may also lead to cross-contamination and therefore infection of a patient. Using an implant that is the wrong size leads to a risk of decreased stress resistance. SPINE INNOVATIONS is not liable for use of another manufacturer's components, or use of its own components for purposes other than those intended. The implants must not be modified or subjected to any treatment.

ADVERSE EFFECTS:

The surgeon is responsible for providing the patient with complete information before the surgical procedure, and especially for informing them of the following potential risks and adverse effects:

- Implant subsidence
- End-plate mismatch
- Pain
- Dysphonia
- Uncomfortable swallowing
- Implant malpositioning (too anterior)
- Retropharyngeal oedema
- Heterotopic ossification
- Allergy
- Horner's syndrome

- Reoperation
- Fusion
- Blood circulation problems
- Degeneration of adjacent levels
- Radiculopathy
- Haematoma
- Disc herniation
- Lower back, leg, hip, or knee problems
- Formation of osteophytes
- Myelopathy
- Infection
- Leakage of cerebrospinal fluid
- Abscess
- Asymptomatic anterolisthesis
- Implant breakage
- Cyst
- Vertigo
- Headaches
- Instability/hypermobility
- Loosening
- Loss of mobility
- Metallosis
- Micro-movement/migration
- Abnormal non-painful sense of touch (tingling, prickling, numbness)
- Paralysis
- Neck tension
- Stenosis
- Trauma during the surgical intervention, such as nerve or spinal cord damage, excess bleeding and/or fractures of vertebral bodies (spinal column bones)
- Bone resorption
- Death

INFORMATION TO GIVE TO THE PATIENT:

The patient must be given:

- The information that allows the device to be identified including the name, batch number, UDI, device model, and name, address, and website of the manufacturer.
- The warnings, precautions, or steps to be taken by the patient or a healthcare professional regarding reciprocal interference from external environmental conditions or sources or medical examinations that can reasonably be predicted.
- All information on the predicted lifespan of the device and the necessary follow-up, which is available in the patient leaflet that can be accessed using the QR code on the implant card issued to the patient.
- As well as all other information intended to guarantee the patient's safe usage of the device, which is available in the patient leaflet that can be accessed using the QR code on the implant card issued to the patient.

Patient instructions are accessible via the QR code on the implant card or at the following website: <https://spine-innovations.com/patient-information>

MRI COMPATIBILITY:

Non-clinical tests have demonstrated that the ESP® range of discs is "MR conditional" in accordance with the definitions of the ASTM + F2503-20 standard. A patient with a device in this range can be scanned safely in an MRI system when the following conditions are met:

- Patient has a single ESP® disc prosthesis implanted
- Patients has no thermoregulation disorder (i.e. no alteration of systemic thermoregulation or decreased local thermoregulation) and
- Patients are under controlled conditions (a doctor and a dedicated trained person who can respond immediately to the physiological stress triggered by heat).
- MRI system with a horizontal tunnel and static magnetic field of 1.5 Tesla or 3 Tesla.
- Spatial gradient of magnetic field below or equal to 19T/m.
- System with B0 ± 1 dB 0 / drl below or equal to 48T²/m.
- Only the whole-body transmission/receipt RF coil is used.
- First level controlled operating mode is used i.e. the specific absorption rate averaged over the whole body (WB-SAR) does not exceed 4 W/kg.
- During non-clinical tests, after 15 minutes of continuous acquisition, the ESP® disc produced a maximum temperature increase of 5.0 ± 1.0 °C at 1.5 T with a WB-SAR measured at 3.50 ± 0.81 W/kg and a maximum temperature increase of 3.5 ± 1.0 °C at 3T with a WB-SAR measured at 3.94 ± 0.88 W/kg.
- With a WB-SAR of 4W/kg, the ESP® disc should produce a maximum temperature increase of 5.7 ± 1.8 °C at 1.5T and 3.6 ± 1.3 °C at 3T.

MRI image quality may be compromised if the area of interest is in the same area as the implant.
It may be necessary to manipulate the acquisition parameters to compensate for this artefact.

IMPORTANT INFORMATION FOR THE SURGEON:

Any serious accident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user or the patient has residence.

PACKAGING AND STERILISATION:

The information on the product label ensures manufacturing traceability and indicates which sterilisation method has been used.

IMPLANTS DELIVERED STERILE:

Verify that the packaging elements (peelable pouches) are perfectly sealed and that all elements are undamaged before using the implants. Do not use a product if its packaging is damaged or if the tamper-proof seal has been broken. The product must not be resterilised by any method whatsoever at the risk of altering its mechanical characteristics.

The sterilisation indicator sticker on the outer packaging should be red for radiation sterilisation. This colour may be altered by poor storage conditions such as heat, humidity, light, etc.

In all cases, an indicator showing the pre-sterilisation colour, yellow, may indicate that the product is not sterile, and, in this case, the product must be returned.

Desterilised instruments cannot be returned.

The expiry date is noted in the product label.

DEVICE DISPOSAL:

In the event of explantation, the implant must be sent to a specialised waste handling service so that it can be disposed of properly with respect to the environment and strict sanitation laws. A defective device that has been removed must be decontaminated and returned to the manufacturer.

DATE OF FIRST CE MARKING:

CP-ESP® DISC PROSTHESIS/CEMENTLESS: 12-2023

The reference text is the French text.

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