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Références.

- AAMI TR12, Designing, testing, and labeling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers.
- ANSI/AAMI ST79, Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- ISO 17665, Sterilization of health care products – Moist heat.
- ISO 17664, Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices.

EN – Instructions for use, cleaning, and sterilisation of SPINE INNOVATIONS instruments

Do not hesitate to consult your supplier's sales department if you need further information. This document is not exhaustive, nor is it a technical surgical manual containing detailed implantation information. Always refer to the documents provided.

General information

This leaflet is applicable to all reusable medical devices manufactured and distributed by SPINE INNOVATIONS and/or its distributors and must be read carefully.

Both new and old instruments must be processed with care in accordance with these instructions before use.

These instructions are not applicable to SPINE INNOVATIONS medical devices supplied sterile, which must not be reused.

Instrument cases do not provide a sterile barrier: sterilisation packaging must be used to keep the device sterile.

- The instructions provided in this leaflet have been validated by SPINE INNOVATIONS. It is the health institution responsibility to ensure that reprocessing is carried out using suitable processing products and equipment and that the staff in charge of reprocessing has received proper training to obtain the expected result. The equipment and processes require regular validation and monitoring.
- Ancillary equipment is supplied non-sterile in transport containers. A transfer sheet (specifying the non-sterile status of the ancillary equipment) is supplied with each piece of ancillary equipment for incoming quality control and prior to sterilisation.
- SPINE INNOVATIONS equipment must be used exclusively for SPINE INNOVATIONS prostheses implantation and explantation.

- Instructions for use**
- A surgical procedure is available from the SPINE INNOVATIONS sales department or its distributors to provide surgeons with additional information on the recommended method for use. It is essential to read this document.
 - It is recommended that this equipment not be used directly after steam sterilisation and to wait for the return to room temperature. Operating suite staff is responsible for checking, prior to the procedure, that SPINE INNOVATIONS instruments are compatible with the operating suite equipment available.
 - Ancillary equipment is designed to operate at room temperature.
 - Some ancillary equipment consists of temporary invasive devices such as drills, reamers, milling cutters, screw taps and probes. These instruments must be used with caution. Their use can damage soft tissues or contact areas.
 - This equipment, tested by SPINE INNOVATIONS, or its distributors, must be handled with care. Any damage to the equipment may give rise to potential malfunction.
 - The equipment must only be used for its intended purpose in surgical procedures.

SPINE INNOVATIONS recommends processing all devices in accordance with the manual cleaning or combined manual/automatic cleaning instructions contained in this leaflet.

- Warnings**
- **The following precautions must be observed** by hospital staff working with contaminated or potentially contaminated medical devices. Pointed or sharp devices must be handled with the utmost caution.
 - **Do not use metal brushes or scouring pads** during manual cleaning procedures. Their use could damage the instrument surface and finish. Use soft nylon bristle brushes and cleaning brushes of suitable sizes on devices to be processed.
 - To facilitate the subsequent steps, **do not allow contaminated devices to dry before reprocessing.**
 - Normal saline solution and cleaning/disinfection agents containing aldehyde, mercury, active chlorine, chloride, bromine, bromide, iodine or iodide are corrosive and **must not** be used.
 - **Do not** use mineral oil or silicone lubricants, as they coat micro-organisms, prevent direct contact of the surface with the steam and are difficult to remove.
 - **Automatic cleaning using only a washer/disinfector is not sufficient for spinal surgery instruments.** Thorough manual cleaning is required before an automatic cleaning procedure.
 - pH-neutral cleaning and enzyme agents are recommended and preferred for cleaning reusable SPINE INNOVATIONS devices.
 - When cleaning aluminium devices, only pH-neutral detergents must be used. Contact with strong alkaline detergents or with solutions containing sodium hydroxide, iodine or chlorine must be avoided, as the aluminium may be affected by chemical corrosion and the device may be damaged. The user must always refer to and comply with the instructions provided by the cleaning product manufacturer.
 - N.B.: Drill bits, pins, rasps and sharp instruments may be treated with alkaline detergents and must be inspected thoroughly after processing to ensure that the sharp edges are not damaged.**
- SPINE INNOVATIONS recommends the use of an ISO 15883 compliant washer/disinfector for instrument cleaning and disinfection.

- Instructions**
- Remove excess bodily fluids and tissues on the instruments with a disposable lint-free pad. Place the devices in a tray of distilled water or cover with a damp cloth.
 - Instruments **must** be cleaned within 30 minutes following use to minimise the risk of drying before cleaning.
 - Used instruments **must** be transported to the procurement department in sealed or covered containers to prevent further risk of contamination.
 - Where applicable, instruments with multiple components must be dismantled for effective cleaning. Take care not to mislay any small screws and components.
 - All cleaning agents must be prepared in accordance with the dilution and temperature recommended by the manufacturer. Softened tap water may be used to prepare cleaning agents. It is important to use the recommended temperatures for optimal cleaning agent performance.
 - N.B.: Cleaning solutions must be prepared fresh.**
 - The instruments **must** be removed from the metal or polymer trays when applying manual and/or automatic cleaning procedures. Lids, cases and instrument trays must be cleaned separately, non-sterile single-use implants are an exception to this rule. Plates and screws may remain in the tray or trolley for reprocessing.
 - **The use of hard water must be avoided.** Softened tap water may be used for the initial rinse. The final rinses must be carried out with purified water to remove deposits on the instruments. One or more of the following processes may be used for water purification: ultrafiltration (UF), reverse osmosis (RO), de-ionisation or equivalent.
 - Steam (moist heat) is the validated sterilisation method for SPINE INNOVATIONS instruments.

- Manual disinfection/cleaning procedure**
- Step 1.** Submerge the instruments completely in an enzyme solution (ALKAZYM 0.5%) and leave to soak for 20 minutes. Use a soft nylon bristle brush to scrub the device gently until any visible stains have been removed. Pay particular attention to rough areas, cavities, contact surfaces, connectors and other areas that are difficult to clean. Move hinged parts to access areas that are difficult to access.
- Step 2.** Remove the device from the enzyme solution and rinse with tap water for at least 3 minutes. Rinse any cavities, holes and any other areas that are difficult to access thoroughly and effectively with a pressurised spray.
- Step 3.** Place the cleaning agents in the ultrasonic cleaner. Submerge the device completely in the cleaning solution and sonicate for at least 10 minutes at 45-50 kHz.
- Step 4.** Rinse the device in purified water for at least 3 minutes or until the blood or visible stains on the device or in the rinsing water have been removed. Rinse any cavities, holes, and any other areas that are difficult to access thoroughly and effectively.
- Step 5.** Repeat the above ultrasonic cleaning and rinsing steps.
- Step 6.** Remove any excess moisture on the instrument with a clean, absorbent, lint-free cloth.

- Manual/automatic disinfection/cleaning procedure**
- Step 1.** Submerge the instruments completely in an enzyme solution (ALKAZYM 0.5%) and leave to soak for 10 minutes. Use a soft nylon bristle brush to scrub the device gently until any visible stains have been removed. Pay particular attention to rough areas, cavities, contact surfaces, connectors and other areas that are difficult to clean. Move hinged parts to access areas that are difficult to access.
- Step 2.** Remove the devices from the enzyme solution and rinse with purified water for at least 3 minutes. Rinse any cavities, holes and any other areas that are difficult to access thoroughly and effectively.
- Step 3.** Place the instruments in a suitable washer/disinfector basket and process via standard instrument washer/disinfector cycle in accordance with the washer/disinfector manufacturer's instructions:
- Wash for 10 minutes at 93°C minimum with a suitable washer/disinfector detergent solution.

- Rinse with demineralised water.
- Dry.

- Inspection before sterilisation**
- The cleaning procedures must be repeated until instruments are "thoroughly cleaned" and "visually clean".
 - Inspect each device carefully to ensure that any visible contamination has been removed. If contamination is observed, repeat the cleaning/disinfection procedure.
 - Check that moving parts (e.g. hinges, connectors, sliding parts, etc.) can move to their full extent.
 - Check that there is no instrument deformation
 - When the instruments have been assembled, check that the devices are assembled correctly with the components in contact.
 - Inspect devices for wear and tear: the device is not suitable for reuse if it looks corroded, discoloured, stained or damaged. Contact your sales representative if the device is defective.

- Sterile packaging**
- Commercially available medical grade steam sterilisation wrapping or pouches may be used to package individual instruments. The wrapping must be prepared according to the DIN EN ISO/ANSI AAMI ISO 11607-1/2 protective double wrapping technique. Trays and cases with lids may also be placed in an approved sterilisation container with a sterilisation sealing lid. Follow the sterilisation container manufacturer's instructions for inserting and replacing sterilisation filters in sterilisation containers.
- Instrument transport and storage cases and trays must be used under the following conditions:**
- If possible, all devices must be disassembled before they are placed in the case.
 - All devices must be arranged to ensure the steam comes into contact with all the instrument surfaces. Instruments must not be stacked or placed in close contact (use the storage holders or silicone mats provided).
 - The user must check that there has been no spillage of the instrument case contents once the devices have been arranged in the case.
 - Only devices manufactured and/or distributed by SPINE INNOVATIONS may be placed in the SPINE INNOVATIONS instrument tray.

- Sterilisation**
- Reusable surgical instruments must always be disinfected according to the specified procedure before sterilisation.
 - The health institution is responsible for in-house instrument reassembly, inspection and packaging procedures following thorough cleaning so as to ensure sterilising steam penetration and suitable drying. Protection measures to be taken against sharp or potentially dangerous instrument edges must also be recommended by the health institution.
 - Steam/moist heat sterilisation is the preferred method recommended for SPINE INNOVATIONS spine surgical instrument sets.
 - The steriliser manufacturer's recommendations must be followed at **all times**. When sterilising several instrument sets in a sterilisation cycle, take care not to exceed the maximum load specified by the manufacturer.
 - The instrument sets must be prepared and packaged correctly in trays and/or cases so that the steam can penetrate and come into direct contact with all the surfaces.
 - See the table below for the cycle sterilisation settings validated by SPINE INNOVATIONS to provide a sterility assurance level (SAL) of 10⁻⁶.
 - Ethylene oxide or plasma sterilisation methods **must not** be used.

Validated steam sterilisation settings:					
	Cycle type	Temperature (°C)	Temperature (°F)	Exposure time	Drying time
A	Prior vacuum	132°C	269.6°F	4 minutes	30 minutes
B	Prior vacuum	134°C	273.2°F	18 minutes	30 minutes
C	Prior vacuum	134°C	273.2°F	3 minutes	30 minutes

Cycles B and C are not recommended in the United States, nor for the inactivation of prions. This procedure was validated by an independent laboratory in accordance with AAMI TIR12.

- Storage**
- When handling packaged cases, take care to avoid damaging the sterile barrier. The health institution must define a storage shelf-life for packaged instruments based on the sterile packaging and the packaging manufacturer's recommendations.

- Responsibilities of the health institution with regard to SPINE INNOVATIONS instruments on loan / Number of potential reuses:**
- SPINE INNOVATIONS does not claim a number of uses of the instruments. If handled and used properly, Spine surgery instruments can be used for several years. However, this service life may be reduced rapidly due to poor handling or unsuitable protection. Instruments that are no longer up to standards due to a long period of use, poor handling or unsuitable care must be returned to SPINE INNOVATIONS. Notify your SPINE INNOVATIONS representative of any problems with the instruments.

- Before being returned to SPINE INNOVATIONS, SPINE INNOVATIONS instruments are required to undergo all the decontamination, cleaning, disinfection, inspection and final sterilisation steps. The decontamination and sterilisation documentation must be provided with the instruments returned to SPINE INNOVATIONS.

IMPORTANT INFORMATION FOR HEALTH CARE ESTABLISHMENT PERSONNEL

Any serious incident related to the device must be reported to the manufacturer and to the competent authority of the Member State of the user's residence.

- References**
- AAMI TIR12, Designing, testing, and labelling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers.
 - ANSI/AAMI ST79, Comprehensive guide to steam sterilisation and sterility assurance in health care facilities.
 - ISO 17665, Sterilisation of health care products – Moist heat.
 - ISO 17664, Sterilisation of medical devices – Information to be provided by the manufacturer for the processing of resterilisable medical devices.

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Referenzen

- AAMI TIR12, Designing, testing, and labeling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers.
- ANSI/AAMI ST79, Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
- ISO 17665, Sterilization of health care products – Moist heat.
- ISO 17664, Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices.

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Bibliografia

- AAMI TIR12. *Designing, testing, and labeling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers.*
- ANSI/AAMI ST79. *Comprehensive guide to steam sterilization and sterility assurance in health care facilities.*
- ISO 17665. Sterilization of health care products – Moist heat.
- ISO 17664. Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices.

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Referenties

- AAMI TIR12, Designing, testing, and labelling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers.
- ANSI/AAMI S779, Comprehensive guide to steam sterilisation and sterility assurance in health care facilities.
- ISO 17665, Sterilisation of health care products – Moist heat.
- ISO 17664, Sterilisation of medical devices – Information to be provided by the manufacturer for the processing of resterilisable medical devices.

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- AAMI TIR12, Designing, testing, and labelling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers.
- ANSI/AAMI ST79, Comprehensive guide to steam sterilisation and sterility assurance in health care facilities.
- ISO 17665, Sterilisation of health care products – Moist heat.
- ISO 17664, Sterilisation of medical devices – Information to be provided by the manufacturer for the processing of resterilisable medical devices.

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Reference

- AAMI TIR12, Designing, testing, and labelling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers.
- ANSI/AAMI ST79, Comprehensive guide to steam sterilisation and sterility assurance in health care facilities.
- ISO 17665, Sterilisation of health care products – Moist heat.
- ISO 17664, Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices.

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Referenser

- AAMI TIR12, Designing, testing, and labelling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers.
- ANSI/AAMI ST79, Comprehensive guide to steam sterilisation and sterility assurance in health care facilities.

- ISO 17665, Sterilisation of health care products – Moist heat.
- ISO 17664, Sterilisation of medical devices – Information to be provided by the manufacturer for the processing of resterilisable medical devices.

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