DISTIMP – IFU-ACB Rev 02 – 02/2025
Cancels and replaces IFU-ACB Rev 01

Latest revision: February 2025
First EC marking date: 2005

Important: The instructions for use are dematerialized and can be accessed via the eIFU web-site link or the QR code provided on the product label. Please ensure you download or print the instructions for use before discarding the label, or keep the label. All electronic versions of the instructions for use can be accessed at: <https://spineWAY.com/en/e-ifu/pros>.

Manufactured by:

CE 2803

1. DEVICE DESCRIPTION AND COMPOSITION

This device is an anterior intersomatic cervical cage. It is intended to reduce the pathologies of the cervical spine by restoring the disc height as well as the physiological bending of the spine and by allowing the fusion of vertebral bodies between them.

This instruction of use applies to the following **ACIFBOX** anterior cervical cages system:

- Cervical cage
- Secure cervical cage
- Screwed cervical cage
- Plates
- Cervical screws
- Product-specific instrumentation (implant holder, trials, square point, screwdriver, etc.)

These components must always be implanted and/or used together. This will ensure that the various component sizes and materials are compatible.

These products are for **single use** and are made available on the market **non-sterile**.

The implants are made of:

- **Cervical cage** - PEEK-OPTIMA® (ASTM F2026) + GOLD (ASTM B562)
- **Secure cervical cage** - PEEK-OPTIMA® (ASTM F2026) + Gold (ASTM B562) + TA6V ELI (EN ISO 5832-3 ASTM F136)
- **Screwed cervical cage** - PEEK-OPTIMA® (ASTM F2026) + Gold (ASTM B562) + TA6V ELI (EN ISO 5832-3 ASTM F136)
- **Plate** – TA6V ELI (EN ISO 5832-3 ASTM F136)
- **Cervical screws** – TA6V ELI (EN ISO 5832-3 ASTM F136)

Use of these components with devices other than those recommended by DISTIMP is prohibited.

The instrumentation is made out of non-implantable stainless steel, TA6V ELI, silicone, HSPF, and PPSU.

2. INDICATIONS FOR USE

Intended population: ACIFBOX cages can be used in skeletally mature individuals.

This medical device is used (and only) to treat cervical disc diseases defined as refractory radiculopathy (irradiant pain) and/or myelopathy (weakness) with herniated disc and/or formation of osteophyte and/or spinal cord compression. The ACIFBOX cervical anterior spinal cages are anatomically designed. They make it possible to restore disc height and physio-logical lordosis. It has graft window intended to be filled with bone substitute, wide enough to achieve a good bone graft.

3. USER

These medical devices must only be implanted and/or used by an HCP who is well-trained in spine surgery. Implantation must be performed using appropriate instrumentation provided by the manufacturer.

4. CONTRAINDICATIONS

The following is a non-exhaustive list of contraindications:

- Acute or chronic, local or systemic infection
- Allergy or intolerance
- Any concurrent disease that could affect implant function, including but not limited to severe osteoporosis, cancer, kidney dialysis, osteopenia, obesity
- Important activity level or impaired mental capacity may be relative contraindications to this surgery.

This medical device is designed, intended and sold only for the uses indicated.

5. POTENTIAL RISKS

- The side effects are the same as those encountered during any cervical spine surgery procedure: infection, pain, hematoma, neuropathy, etc.
- Intraoperative neurological complications which may require temporary or permanent removal of the fixation hardware.
- Secondary neurological complications which may require a revision procedure for partial or complete, temporary or permanent removal of the fixation hardware.
- Allergic reaction to implant materials
- Pain, discomfort due to the presence of the device
- Decrease in bone sensitivity due to stress shielding
- Components loosening or disassembly, requiring an additional surgical procedure
- Non-union (pseudarthrosis) or delayed union
- Bending or fracture due to mechanical wear over time.
- In case of looseness between the components, a disassembly may appear and may require additional intervention.
- Due to the anterior cervical surgical approach, the following effects may occur hoarseness or difficulty swallowing, disease of the adjacent segment, nerve damage.

Warning: Patients receiving ACIFBOX cage should be advised that implant longevity may be affected by their weight, age and activity level.

6. WARNINGS

- The patient has to be informed of the risks of the surgical intervention by the medical staff.

The patient should follow the advice and recommendations from the surgeon (radiological examinations) during the postoperative phase.

- It is recommended to limit physical activity after the surgery to allow bone fusion. If not, the implant may break or otherwise be damaged necessitating revision surgery. The implant should not be exposed to extreme movements of mechanical vibrations for example.
- All modifications (aspect, pain...) at the implant site must be reported to the medical practitioner. The medical practitioner should also be informed of all type of incidents like a fall for instance, even when there are no visible signs at the implant site.
- Never reuse a component of the ACIFBOX system that has previously been implanted. Reuse of the device is prohibited because of the chemical, biological (allergy, toxicity, contamination, infection) and mechanical (deterioration, implant wear, etc.) risks.
- The use of nonsteroidal anti/inflammatory drugs (NSAIDs) may be limited due to interferences and delays in the bone healing process.
- **PATIENT SELECTION** will greatly affect the results.
 - Patients suffering from obesity, malnutrition, and/or poor bone quality and patients who smoke, or abuse alcohol are poor candidates for spinal fusion. Patients who smoke should be advised of the consequences of the fact that an increased incidence of non-union has been reported when smoking
 - Obese patient presents a significant increased morbidity profile.
 - Patients with previous spinal surgery at the level(s) to be treated may have different clinical outcomes compared to those without previous surgery.
 - An active, debilitated, or demented patient who cannot properly follow recommendations may be particularly at risk during postoperative rehabilitation.

7. PRECAUTIONS

1. In case of end plate damage, the implantation of an additional stabilization system might be recommended.
2. To avoid an encroachment of bone fusion mass into the spinal canal, in case of narrow canal, the use of a stand-alone cage group should be limited.
3. Do not use any damaged or suspect device, even if it is new.
4. **ADEQUATELY INSTRUCT THE PATIENT.** Postoperative care and the patient's ability and willingness to follow instructions are among the most important aspects of successful bone healing. Inform the patient about the implant limitations, the expected bone healing achievement time frame and the need for avoiding any other activity that would compromise or delay the healing process, such as smoking, consuming alcohol and inadequate physical activities, especially lifting and twisting motions. Patient must be informed that an implant is not as strong as normal healthy bone or disc and could loosen, subside and/or break in case of over mechanical solicitations, especially in the absence of complete bone healing. Implants displaced or damaged by improper activities may migrate and damage the neurological structures or blood vessels.
5. ACIFBOX spinal implants have a finite useful life. They are not intended or expected to be the only mechanism for spinal support. Without bony fusion, the devices are not expected to support the spine forces indefinitely and will fail in any of several modes. Even if solid bone fusion occurs, implant components may nevertheless subside, break, or loosen. Therefore, the patient must be made aware that implant components may subside, break, or loosen even though restrictions in activity are followed.
6. **IMPLANT REMOVAL.** The implant is intended to stay in place after performance achievement. However, the surgeon may remove these implants if bone fusion doesn't occur. The possibility of a second surgical procedure and the risks associated with this second surgical procedure must be discussed with the patient and weighed versus benefits. If the implants do break, or in any other malfunction, the decision to remove them must be made by the physician who must consider the condition of the patient and the associated risks.

8. RISKS OF INTERFERENCE IN MEDICAL IMAGING EXAMINATIONS (MRI/CT SCAN)

Security in magnetic resonance environment (MRI):

The ACIFBOX System is MR Conditional. A person with these devices may be safely scanned under the following conditions:

Static Magnetic Field Strength (Bo)	1.5 or 3.0T
Maximum Spatial Field Gradient	7.0 T/m (700 gauss/cm)
RF Polarization	Circularly Polarized (CP)
RF Transmit Coil Type	Body coil
RF Operating Mode	Normal Operating Mode
Maximum Whole-Body SAR	2 W/kg
Scan Duration	20 minutes of scanning (per pulse sequence)

Heating:

The maximum temperature rise measured is 3.8 °C after 15 minutes of continuous scanning in a 1.5T scanner with a MR sequence with a SAR of 2 W/kg. Scanning in a 3T scanner generates lower heating.

Based on the temperature versus time data, a cooling time of 10 minutes was determined to be sufficient for the temperature of each device, after the scan sequence was stopped, to reach a steady state.

Disclaimer:

The MR safety has been evaluated in the above-mentioned conditions. The MR safety in other conditions has not been evaluated and thus is unknown. MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the implant/device.

9. DEVICES SUPPLIED NON STERILE – CLEANING, DECONTAMINATION AND STERILIZATION

All instruments and implants provided non-sterile must first be thoroughly cleaned using established hospital methods before sterilization and introduction into a sterile surgical field, per the directions below.

Additionally, all instruments that have been previously taken into a sterile surgical field must be immediately decontaminated and thoroughly cleaned using established hospital methods before sterilization and reintroduction into the sterile surgical field.

Repeated processing cycles that include automated cleaning and steam sterilization have minimal effects on the devices. For users who do not use the recommended methods, we recommend that they validate the methods they use by appropriate laboratory techniques.

ACIFBOX - Cervical cage - Instruction for use

Note1: Certain cleaning solutions such as those containing sodium hypochlorite or formalin may damage devices, particularly instruments.

Note2: In the case of high-risk patients being suspected of contact with Unconventional Transmissible Agents or prions (e.g., Creutzfeldt- Jacob disease), decontamination should be performed according to the recommendations of the World Health Organization. In this case, DISTIMP instruments can be decontaminated with sodium hydroxide.

9.1 Decontamination at the point of use:

Where possible, instruments shall be disassembled.

This operation is MANDATORY for dirty instruments to lower the microorganism population and make easier the further cleaning.

This operation shall be done with gloves 30 minutes after use to limit the possibility of soil drying.

- Soak instruments in a neutral pre-disinfection bath in order to avoid corrosion and a temperature lower than 30°C to prevent microorganism fixation.
 - Articulated instruments (ex: clamp, scissors,) shall be opened.
 - Instruments with several parts shall be dismantled.
 - During handling, a particular attention shall be paid on sharp instruments.
 - Soaking time shall comply with soaking manufacturer indications.
 - Instruments must be rinsed out with running water (temperature lower than 30°C).
- NOTE: the pre-disinfection bath must be replaced after each operation.

9.2 Cleaning:

Automated cleaning process has been validated according to ISO 17664-1 and AAMI TIR 30.

The automated cleaning cycle to be used is as follows:

1. Pre-cleaning - temperature < 45°C - duration 2min minimum
2. Cleaning - temperature ≈ 55°C - duration 5min minimum
3. Neutralizing - approximate duration 2min
4. Rinsing with cold tap water - approximate duration 2min
5. Thermal Disinfection with deionized water - temperature ≈ 90°C - duration 5min minimum
6. Drying of the devices with a sterile gauze pad. Clean compressed air may be used to dry the lumens.
7. Perform a visual inspection on the devices and verify that they are clean, dry, and in proper working order prior to sterilization.

The automated cleaner need to be qualified according to the requirements defined in the ISO15883 standards.

Cleaning solution (detergent) has to be selected by the hospital to be compatible with the use in automated cleaner/disinfector according to the ISO 15883 standard and on surgical instruments made out of metallic and/or plastic parts.

Doses, temperatures, and durations may be adapted according to solution manufacturer instructions.

Because of the use of some cleaning solutions, white traces may appear after drying. A discoloration has no adverse effect on titanium alloy devices.

9.3 Inspection and functional test:

Instruments:

No maximal number of cycles is recommended but, between each use, all instruments need to be:

- Visually inspected for trace of organic residues corrosion (rust, pitting), damages such as scratches and notches, debris, discoloration, residue, flaking, wear, cracks and unreadable laser markings.
- Test for proper functioning of all joints, moving pieces and mechanisms and test the functionality of assemblies and connections. When needed, proceed to maintenance before sterilization.

Where appropriate, instruments are designed to be reassembled the identical way they were disassembled.

Any damaged instruments (for example with corrosion, excessive scratches, notches, residue, debris or unreadable laser markings) should be destroyed and replaced by a new instrument. In case of unexpected damage, please return the device to DISTIMP as a complaint.

Before sterilization check that the instruments are all there and are in good working order.

Implants and instruments:

For sterilization and storage, the devices can be settled in metallic containers and arranged in kits following carefully the DISTIMP recommendations to assure sterilization efficacy. Then, containers can be sterilized with the devices.

9.4 Sterilization:

All DISTIMP instruments and implants that are provided non-sterile and must be steam sterilized by hospitals and clinics according to one of the following validated methods.

Method	Cycle	Temperature	Exposure time	Drying time
Steam	Vacuum	134°C	18 min. minimum	-
Steam	Pre-vacuum	270°F (132°C)	4 min.	20 min.
Steam	Pre-vacuum	134°C	3 min.	20 min.

Note for France: Use the first cycle parameters according to the French Circular DGS/R13/2011/449 of 1 December 2011.

These methods are validated through ANSI/ AAAMI ST79 to meet the requirements of a SAL of 10⁻⁶.

If after having followed this sterilization method there is still water in the sterilization containers or on/inside the device, the device must be dried, and sterilization repeated.

9.5 Maintenance:

The instruments joints, hinges, moving pieces and mechanisms should be regularly lubricated with a lubricant meant for medical grade surgical instruments, to reduce friction and wear. Complying with the lubricant supplier's instructions is highly recommended. The lubricant used must be compatible with steam sterilization.

10. STORAGE AND DISPOSAL

- Store devices in a manner that provides protection from dust, light, insects, vermin and extreme temperature and humidity.
- Avoid contact with chemical products or corrosive vapors.
- Do not store the instrument close to products that may have a corrosive action (chlorine).
- Unpacked devices can be stored in the dedicated container.
- It is IMPERATIVE to distinguish STERILE and NON-STERILE instruments.

The disposal of products shall be done according to the current procedures of the health care center to prevent any risk of cross-contamination.

11. WARRANTY

All warranty rights are lost if repairs or modifications are carried out by an unauthorized service center.

The manufacturer does not take responsibility for any effects on safety, reliability or performance of the product if the product is not used in conformity with the instructions for use.

Technical alterations reserved.

12. COMPLAINTS/ INCIDENT AND INFORMATION

Any serious incident that occurs in relation to the device must be reported to the manufacturer as a complaint, and for Europe to the national competent authority of the patient/surgeons.

Any other grounds for dissatisfaction relating to the quality of the product should be notified to DISTIMP or its representative.

For all complaints, please give the name and reference along with the batch number of the component(s), your name and address and an exhaustive description of the event to help DISTIMP understand the cause of the complaint. When possible, please send the device to DISTIMP for the investigation.

All current information regarding ACIFBOX products (instruction for use, surgical technique manual, check-list) is available on the website <https://www.distimp.com>.

The instructions for use are available in paper form upon request at no additional cost within 7 calendar days.

DISTIMP SAS - 7 allée du Moulin Berger
69130 ECULLY, FRANCE
Tel : +33 (0)4 72 77 01 52
www.distimp.com
quality@distimp.com
info@distimp.com

SYMBOL KEY LEGEND FROM ISO 15223-1

Caution: symbol key is for reference only - some symbols listed may not apply. Please see main product label affixed to product.

MANUFACTURER

NON-STERILE

DO NOT REUSE

KEEP DRY

CATALOG NUMBER

KEEP AWAY FROM SUNLIGHT

CAUTION: CONSULT
ACCOMPANYING DOCUMENTS

UNIQUE DEVICE IDENTIFIER

BATCH NUMBER

MR CONDITIONAL

CONSULT INSTRUCTIONS
FOR USE

MEDICAL DEVICE

MANUFACTURING DATE

DO NOT USE IF PACKAGE IS
DAMAGED AND CONSULT
INSTRUCTIONS FOR USE

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