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the device, the instrument must be dried and sterilization repeated.

- Do not use instruments if the functional parts are damaged.

#### 7) PROCESSING INSTRUCTIONS

Instruments can go through this cycle a number of times (read carefully the inspection section).

For users who do not use the recommended method, we recommend that they validate the methods they use by appropriate laboratory techniques.

Instruments shall be decontaminated immediately after each use. They shall be thoroughly cleaned before sterilization and reintroduction into the sterile surgical field.

#### DECONTAMINATION AT THE POINT OF USE

Where possible, instruments shall be disassembled.

This operation is MANDATORY for dirty instruments to lower the microorganism population and make easier the further cleaning. This operation shall be done with gloves 30 minutes after use to limit the possibility of soil drying.

- Soak instruments in a neutral pre-disinfection bath in order to avoid corrosion and a temperature lower than 30°C to prevent microorganisms fixation.
- Articulated instruments (ex: clamp, scissors,...) shall be opened.
- Instruments with several parts shall be dismantled.
- During handling, a particular attention shall be paid on sharp instruments.
- Soaking time shall comply with manufacturer indications.
- Instruments must be rinsed out with running water (temperature lower than 30°C).

NOTE: The pre-disinfection bath must be replaced after each operation.

- Unpacked instruments can be stored in the dedicated container.
- It is IMPERATIVE to distinguish STERILE and NON-STERILE instruments.

The disposal of products shall be done according to the current procedures of the health care center to prevent any risk of cross-contamination.

**10) COMPLAINTS/ INCIDENT AND INFORMATION**
Any serious incident that occurs in relation to the device must be reported to the manufacturer as a complaint, and for Europe to the national competent authority of the patient/surgeons.

Any other grounds for dissatisfaction relating to the quality of the product should be notified to Spineway or its representative.

For all complaints, please give the name and reference along with the batch number of the component(s), your name and address and an exhaustive description of the event to help Spineway understand the cause of the complaint. Please send the device to Spineway for the investigation.

All current information regarding Spineway products (instruction for use, surgical technique manual, check-list) is available on the website https://www.spineway.support and on the MySpineway App (from the App Store or Play Store), section Products.

For further information, complaints, or to obtain Surgical Technique Manuals please contact Spineway at the following address or your representative.

**CLEANING**
Automated cleaning process has been validated according to ISO 17664-1 and AAMI TIR 30. The automated cleaning cycle to be used is as follows:

1. Pre-cleaning - temperature < 45°C - duration 2min minimum

2. Cleaning - temperature = 55°C - duration 5min minimum

3. Neutralizing - approximate duration 2min

4. Rinsing with cold tap water - approximate duration 2min

5. Thermal Disinfection with deionized water - temperature = 90°C - duration 5min minimum

6. Drying of the instruments with a sterile gauze pad. Clean compressed air may be used to dry the lumens.

7. Perform a visual inspection on the instruments and verify that they are clean, dry, and in proper working order prior to sterilization.

The automated cleaner needs to be qualified according to the requirements defined in the ISO15883 standards.

Cleaning solution (detergent) has to be selected by the hospital to be compatible with the use in automated cleaner/disinfector according to the ISO 15883 standard and on surgical instruments made out of metallic and/or plastic parts. Using an alkaline reagent is recommended, such as Neodisher.

Doses, temperatures and durations may be adapted according to solution manufacturer instructions. Because of the use of some cleaning solutions, white traces may appear after drying. For the black aluminum dilator previously marketed, do not use an alkaline reagent.

A disoloration has no adverse effect on titanium alloy devices.

#### INSPECTION AND FUNCTIONAL TEST

No maximal number of cycles is recommended but, between each use, all instruments to be used:

- visually inspected for trace of organic residues, corrosion (rust, pitting), damages such as scratches and notches, debris, discoloration, residue, flaking, wear, cracks and unreadable laser markings

-test for proper functioning of all joints, moving pieces and mechanisms and test the functionality of assemblies and connections. When needed, proceed to maintenance before sterilization.

Additional test for the Guiding wires: the straightness needs to be visually and functionally inspected (the guiding wire shall pass through the cannula of the screw).

Cannulated taps should be checked as some material could obstruct the canula after drilling. A guiding wire should be inserted into the drill canula to eliminate any some material and avoid obstruction.

Where appropriate, instruments are designed to be reassembled the identical way they were disassembled. Any damaged instruments (for example with corrosion, excessive scratches, notches, residue, debris or unreadable laser markings) should be destroyed and replaced by a new instrument.

For sterilization and storage, the instruments can be settled in metallic containers and arranged in kits following carefully the Spineway recommendations to assure sterilization efficacy. Then, containers can be sterilized with the instruments.

##### STERILIZATION

All Spineway instruments are provided non-sterile and must be steam sterilized by hospitals and clinics according to one of the following validated methods:

Method	Cycle	Temperature	Exposure time	Drying
Steam	Vacuum	134°C	18 min. minimum	-
Steam	Pre-vacuum	220° (132°C)	4 min.	20 min.
Steam	Pre-vacuum	134°C	3 min.	20 min.

Note for France: Use the first cycle parameters according to the French Circular DGS/R3/2011/449 of 1 December 2011.

Note for the USA: The user is instructed to use FDA cleared wraps for the sterilization process.

These methods are validated through ANSI/ AAMI ST79 to meet the requirements of a SAL of 10<sup>-6</sup>.

If after having followed this sterilization method there is still water in the sterilization containers or on/inside

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**8) MAINTENANCE**
The instruments joints, hinges, moving pieces and mechanisms should be regularly lubricated with a lubricant meant for medical grade surgical instruments, to reduce friction and wear.
Complying with the lubricant supplier's instructions is highly recommended. The lubricant used must be compatible with steam sterilization.

**9) STORAGE AND DISPOSAL**
Store devices in a manner that provides protection from dust, light, insects, vermin and extremes of temperature and humidity.

- Avoid contact with chemical products or corrosive vapors.
- Do not store the instrument close to products that may have a corrosive action (chlorine).

Unpacked instruments can be stored in the dedicated container.

- It is IMPERATIVE to distinguish STERILE and NON-STERILE instruments.

The disposal of products shall be done according to the current procedures of the health care center to prevent any risk of cross-contamination.

**10) COMPLAINTS/ INCIDENT AND INFORMATION**
Any serious incident that occurs in relation to the device must be reported to the manufacturer as a complaint, and for Europe to the national competent authority of the patient/surgeons.

Any other grounds for dissatisfaction relating to the quality of the product should be notified to Spineway or its representative.

For all complaints, please give the name and reference along with the batch number of the component(s), your name and address and an exhaustive description of the event to help Spineway understand the cause of the complaint. Please send the device to Spineway for the investigation.

All current information regarding Spineway products (instruction for use, surgical technique manual, check-list) is available on the website https://www.spineway.support and on the MySpineway App (from the App Store or Play Store), section Products.

For further information, complaints, or to obtain Surgical Technique Manuals please contact Spineway at the following address or your representative.

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