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For users who do not use the recommended method, we recommend that they validate the methods they use by appropriate laboratory techniques.

Instruments shall be decontaminated immediately after each use. They shall be be thoroughly cleaned before sterilization and reintroduction into the sterile surgical field.

DECONTAMINATION AT THE POINT OF USE

Where possible, instruments shall be disassembled.

This operation is MANDATORY for dirty instruments to lower the microorganism population and make easier the further cleaning.

This operation shall be done with gloves 30 minutes after use to limit the possibility of soil drying.
• Soak instruments in a neutral pre-disinfection bath in order to avoid corrosion and a temperature lower than 30°C to prevent microorganism fixation.

• Articulated instruments (ex: clamp, scissors,...) shall be opened.
• Instruments with several parts shall be dismantled.
• During handling, a particular attention shall be paid on sharp instruments.
• Soaking time shall comply with manufacturer indications.
• Instruments must be rinsed out with running water (temperature lower than 30°C).
NOTE: the pre-disinfection bath must be replaced after each operation.

Note: In the case of high-risk patients being suspected of contact with Unconventional Transmissible Agents or prions (e.g. Creutzfeldt-Jacob disease), decontamination should be performed according to the recommendations of the World Health Organization. In this case, Distimp products can be decontaminated with sodium hydroxide.

CLEANING

Automated cleaning process has been validated according to ISO 17864-1 and AAMI TR30.

The automated cleaning cycle to be used is as follows:

- Pre-cleaning - temperature < 45°C - duration 2min minimum
- Cleaning - temperature = 55°C - duration 5min minimum
- Neutralizing - approximate duration 2min
- Rinsing with cold tap water - approximate duration 2min
- Thermal disinfection with deionized water - temperature = 90°C - duration 5min minimum
- Drying of the instruments with a sterile gauze pad. Clean compressed air may be used to dry the lumen.
- Perform a visual inspection on the instruments and verify that they are clean, dry, and in proper working order prior to sterilization.

The automated cleaner needs to be qualified according to the requirements defined in the ISO15883 standards.

Cleaning solution (detergent) has to be selected by the hospital to be compatible with the use in automated cleaner/ disinfecter according to the ISO 15883 standard and on surgical instruments made out of metallic and/or plastic parts. Using an alkaline reagent is recommended, such as Neodisher.

Doses, temperatures and durations may be adapted according to solution manufacture instructions. Because of the use of some cleaning solutions, white traces may appear after drying. For the black aluminium drain previously marketed, do not use an alkaline reagent.

A discoloration has no adverse effect on titanium alloy devices.

INSPECTION AND FUNCTIONAL TEST

No maximal number of cycles is recommended but, between each use, all instruments need to be:
visually inspected to trace of organic residues, corrosion (rust, pitting), damages such as scratches and notches, debris, discoloration, residue, fishing, wear, cracks and unreadable laser markings.

-test for proper functioning of all joints, moving pieces and mechanisms and test the functionality of assemblies and connections. When needed, proceed to maintenance before sterilization.

Additional test for the Guiding wires: the straightness needs to be visually and functionally inspected (the guiding wire shall pass through the cannula of the screw).

Cannulated taps should be checked as bone material could obstruct the canula after drilling. A guiding wire should be inserted into the drill canula to eliminate any bone material and avoid obstruction.

Where appropriate, instruments are designed to be reassembled the identical way they were disassembled.

Any damaged instruments, which may compromise the function of the instrument (for example with corrosion, excessive scratches, instrument's shape deformation, notches, residue, debris or unreadable laser markings) should be discarded and replaced by a new instrument.

For sterilization and storage, the instruments can be settled in metallic containers and arranged in kits following carefully the Distimp recommendations to assure sterilization efficacy. Then, containers can be sterilized with the instruments.

STERILIZATION

As distimp instruments are provided non-sterile and must be steam sterilized by hospitals and clinics according to one of the following validated methods.

Méthode	Cycle	Température	Exposition time	Drying time
Steam	Vacuum	134°C	18 min. minimum	-
Steam	Pre-vacuum	220°F (102°C)	4 min.	20 min.
Steam	Pre-vacuum	134°C	3 min.	20 min.

Note for France: Use first cycle parameters according to the French Cirular DGS/R3/2011/449 of 1 December 2011.

Note for the USA: The user is instructed to use FDA cleared wars for the sterilization process.

These methods are validated through ANSI/ AAMI ST79 to meet the requirements of a SAL of 10⁻⁶.

If after having followed this sterilization method there is still water in the sterilization containers or on/inside the device, the device must be dried and sterilization repeated.

8) MAINTENANCE

The instruments joints, hinges, moving pieces and mechanisms should be regularly lubricated with a lubricant meant for medical grade surgical instruments, to reduce friction and wear.
Complying with the lubricant supplier's instructions is highly recommended. The lubricant used must be compatible with steam sterilization.

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9) STORAGE AND DISPOSAL
Store devices in a manner that provides protection from dust, light, insects, vermin and extremes of temperature and humidity.
-Avoid contact with chemical products or corrosive vapors.
-Do not store the instrument close to products that may have a corrosive action (chlorine).

Unpacked instruments can be stored in the dedicated

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