IFU-VLS Rev 01  
Latest revision: July 2024

**VEOS - Spine fixation system - Instruction for use**

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**CE 2803**

DISTIMP – IFU-VLS Rev 01 – 07/2024  
Cancels and replaces IFU-VLS Rev 00

Latest revision: July 2024  
First EC marking date: 2008

Manufactured by:

CE 2803

## 1. DEVICE DESCRIPTION AND COMPOSITION

This device is a posterior osteosynthesis system. It is intended for the reduction of pathologies of the thoraco-lombo-sacral spine by restoring the disc height as well as the physiological bending of the spine.

This instruction of Use applies to the following **VEOS** spine fixation system components:

- Monoaxial, polyaxial and polyaxial spondylo screws
- Bended and straight titanium rod
- Bended PEEK rod
- Connectors
- Transverse connector
- Nut
- Product-specific instrumentation (forceps, rod bender, rod pusher, screwdriver, etc.)

Before the surgery, read and refer to the surgical technique manuals to support you in selecting suitable implants, instruments and for a safe combination of these devices. Surgical technique manuals provide also instructions to use devices safely and to verify whether implants are properly installed. A good understanding of the surgical technique manuals is essential to also prevent misuse that may cause injury to the patient or users.

**CORRECT IMPLANT SELECTION IS VITAL.** Selecting the proper implant size, shape, and design increases the potential for satisfactory fixation.

These products are for **single use** and are made available on the market **non-sterile**.

The implants are made of TA6V titanium alloy (ISO 5832-3, ASTM F136) or PEEK-OPTIMA® LT1-6BA (ASTM F2026).

Use of these components with devices other than those recommended by DISTIMP is prohibited. The instrumentation is made out of non-implantable stainless steel, silicone, nitinol, HSPP and PPSU.

## 2. INDICATIONS FOR USE

Intended population: VEOS system can be used in skeletally mature individuals (adults) and in adolescent population for specific indications such as scoliosis.

This medical device is intended for the posterior fixation of the thoracic, lumbar and sacral spine. It is indicated for the treatment of deformities (any etiology), trauma, tumors, and the degenerative spine (spondylolisthesis, degenerative disc disease, spinal fractures, spinal stenosis, non-union). VEOS can also be used in revision surgery in case of failure of a previous procedure.

## 3. USER

These medical devices must only be implanted and/or used by an HCP who is well-trained in spine surgery. Implantation must be performed using appropriate instrumentation provided by the manufacturer.

## 4. CONTRAINDICATIONS

The following is a non-exhaustive list of contraindications:

- Acute or chronic, local or systemic infection
- Allergy or intolerance
- Any concurrent disease that could affect implant function, including but not limited to severe osteoporosis, cancer, kidney dialysis, osteopenia, obesity
- Important activity level or impaired mental capacity may be relative contraindications to this surgery.

This medical device is designed, intended and sold only for the uses indicated.

## 5. POTENTIAL RISKS

- The side effects are the same as those encountered during any spine surgery procedure: infection, pain, hematoma, neuropathy, etc.
- Intraoperative neurological complications which may require temporary or permanent removal of the fixation hardware.
- Secondary neurological complications which may require a revision procedure for partial or complete, temporary or permanent removal of the fixation hardware.
- Allergic reaction to implant materials
- Pain, discomfort due to the presence of the device
- Decrease in bone sensitivity due to stress shielding
- Components loosening or disassembly, requiring an additional surgical procedure
- Non-union (pseudarthrosis) or delayed union
- Bending or fracture due to mechanical wear over time.
- Disease of the adjacent segment.

**Warning:** Patients receiving the VEOS system should be advised that implant longevity may be affected by their weight, age and activity level, and that premature or inappropriate physical activity could also reduce its longevity.

## 6. WARNINGS

- Post-operative instructions, depending on the scope of the instrumentation, its level and the pathology involved, must be transmitted to the patient in order to be taken into account by the latter.
- A second procedure may be needed for anterior spine fixation after the posterior fixation: either planned beforehand as a part of a two-stage procedure or because of the patient's anatomy and/or pathology.
- Never reuse a component of the VEOS system that has previously been implanted. Reuse of the device is prohibited because of the chemical, biological (allergy, toxicity, contamination, infection) and mechanical (deterioration, implant wear, etc.) risks.

## 7. PRECAUTIONS OF USE

1. **SURGICAL IMPLANTS MUST NEVER BE REUSED.** Even though the device appears undamaged, it may have small defects and internal stress patterns which may lead to early

breakage. In the event of a revision, organize the surgery in order to systematically replace all the existing components manipulated during surgery.

2. Do not use any damaged or suspect device, even if it is new.
3. **CORRECT IMPLANT HANDLING IS VITAL.** Only contour metal implants with proper equipment. Avoid any notching, scratching or reverse bending of the devices when contouring. Alterations will produce defects in surface finish and internal stresses that may become the focal point for eventual breakage. Do not use the implant if damage is suspected. Bending of screws will significantly decrease the fatigue life and may cause failure. Excessive torque applied to the screws when seating the rod may cause bone failure resulting in stripped threads and/or compromised screw purchase.
4. **BENDING THE CONSTRUCT.** Titanium alloy components should never be bent sharply or reverse bent. If a construct is over-contoured contour a new construct correctly rather than reverse bending the over-contoured construct.
5. Patient selection: the following factors can be extremely important to the success of the procedure:
  - **The patient's weight.** An overweight patient can produce loads on the device that can lead to failure of the appliance and the operation. This profile of patient increases the incidence and/or severity of morbidity in the postoperative course.
  - **Senility, mental illness, alcoholism, or drug abuse.** These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions in the device use, leading to implant failure or other complications.
  - **Certain degenerative diseases.** In some cases, degenerative disease may be so advanced that it may substantially decrease the device's expected performance. For such cases, orthopedic non-surgical treatment (i.e: braces, halo-traction) can only be considered as temporary solution.
  - **Foreign body sensitivity.** No pre-operative test can completely exclude the possibility of sensitivity or allergic reaction. Patients can develop sensitivity or allergy after implants have been in the body for a period of time.
  - **Smoking.** Patients who smoke have been observed to experience delayed healing and higher rates of pseudoarthrosis following surgical procedures where bony fusion is required. Additionally, smoking has been shown to be a risk for complication after surgery like pulmonary and cardiovascular complications.
  - **Patients with previous spinal surgery** at the level(s) to be treated may have different clinical outcomes compared to those without previous surgery for similar indication.
  - **Diabetes mellitus** is associated with increased risk of surgical site infection.
6. **IMPLANT REMOVAL AFTER HEALING.** The surgeon may remove these implants after bone fusion occurs. The possibility of a second surgical procedure and the risks associated with this second surgical procedure must be discussed with the patient and weighed versus benefits. If the implants do break, or in any other malfunction, the decision to remove them must be made by the physician who must consider the condition of the patient and the associated risks. If the device is not removed after the completion of its intended use, any of the following complications may occur: (1) Corrosion, with localized tissue reaction or pain; (2) Implant migration resulting in injury; (3) Risk of additional injury from postoperative trauma; (4) Bending, loosening, and/or breakage, which could make removal impractical or difficult; (5) Pain, discomfort, or abnormal sensations due to device presence; (6) Possible increased risk of infection; (7) Bone loss due to stress shielding.
7. **ADEQUATELY INSTRUCT THE PATIENT.** Postoperative care and the patient's ability and willingness to follow instructions are among the most important aspects of successful bone healing. Inform the patient about the implant limitations, the expected bone healing achievement time frame and the need for avoiding any other activity that would compromise or delay the healing process, such as smoking, consuming alcohol and inadequate physical activities, especially lifting and twisting motions. Patient must be informed that a metallic implant is not as strong as normal healthy bone and could loosen, bend and/or break in case of over mechanical solicitations, especially in the absence of complete bone healing. Implants displaced or damaged by improper activities may migrate and damage the neurological structures or blood vessels. An active, debilitated, or demented patient who cannot properly use weight-supporting devices may be particularly at risk during postoperative rehabilitation.
8. **PHYSICIAN POSTOPERATIVE MOBILIZATION:** Until X-rays confirm fusion mass maturation, external immobilization (such as bracing or casting) is recommended. Instructions to the patient to reduce stress on the implants are an equally important part of the attempt to avoid the occurrence of clinical problems that may accompany fixation failure.
9. Use of VEOS system components with devices other than those recommended by DISTIMP is prohibited.
10. VEOS implants have a finite useful life. They are not intended or expected to be the only mechanism for spinal support. Without bony fusion, the devices are not expected to support the spine forces indefinitely and will fail in any of several modes. Even if solid bone fusion occurs, implant components may nevertheless bend, break, or loosen. Therefore, the patient must be made aware that implant components may bend, break, or loosen even though restrictions in activity are followed.

## 8. RISKS OF INTERFERENCE IN MEDICAL IMAGING EXAMINATIONS (MRI / CT SCAN)

Security in magnetic resonance environment (MRI):

The VEOS Spinal Fixation System is MR Conditional. A person with these devices may be safely scanned under the following conditions:

|                                     |                                             |
|-------------------------------------|---------------------------------------------|
| Static Magnetic Field Strength (Bo) | 1.5 or 3.0T                                 |
| Maximum Spatial Field Gradient      | 7.0 T/m (700 gauss/cm)                      |
| RF Polarization                     | Circularly Polarized (CP)                   |
| RF Transmit Coil Type               | Body coil                                   |
| RF Operating Mode                   | Normal Operating Mode                       |
| Maximum Whole-Body SAR              | 2 W/kg                                      |
| Scan Duration                       | 20 minutes of scanning (per pulse sequence) |

### Heating:

The maximum temperature rise measured is 3.8 °C after 15 minutes of continuous scanning in a 1.5T scanner with a MR sequence with a SAR of 2 W/kg. Scanning in a 3T scanner generates lower heating.

Based on the temperature versus time data, a cooling time of 10 minutes was determined to be sufficient for the temperature of each device, after the scan sequence was stopped, to reach a steady state.

### Disclaimer:

The MR safety has been evaluated in the above-mentioned conditions. The MR safety in other conditions has not been evaluated and thus is unknown.

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the implant/device.

# VEOS - Spine fixation system - Instruction for use

## 9. DEVICES SUPPLIED NON STERILE – CLEANING, DECONTAMINATION AND STERILIZATION

All instruments and implants provided non-sterile must first be thoroughly cleaned using established hospital methods before sterilization and introduction into a sterile surgical field, per the directions below.

Additionally, all instruments that have been previously taken into a sterile surgical field must be immediately decontaminated and thoroughly cleaned using established hospital methods before sterilization and reintroduction into the sterile surgical field.

Repeated processing cycles that include automated cleaning and steam sterilization have minimal effects on the devices. For users who do not use the recommended methods, we recommend that they validate the methods they use by appropriate laboratory techniques.

Note1: Certain cleaning solutions such as those containing sodium hypochlorite or formalin may damage devices, particularly instruments.

Note2: In the case of high-risk patients being suspected of contact with Unconventional Transmissible Agents or prions (e.g. Creutzfeldt-Jacob disease), decontamination should be performed according to the recommendations of the World Health Organization. In this case, DISTIMP instruments can be decontaminated with sodium hydroxide.

### 9.1 Decontamination at the point of use:

Where possible, instruments shall be disassembled.

This operation is MANDATORY for dirty instruments to lower the microorganism population and make easier the further cleaning.

This operation shall be done with gloves 30 minutes after use to limit the possibility of soil drying.

- Soak instruments in a neutral pre-disinfection bath in order to avoid corrosion and a temperature lower than 30°C to prevent microorganism fixation.
- Articulated instruments (ex: clamp, scissors,...) shall be opened.
- Instruments with several parts shall be dismantled.
- During handling, a particular attention shall be paid on sharp instruments.
- Soaking time shall comply with soaking manufacturer indications.
- Instruments must be rinsed out with running water (temperature lower than 30°C).

NOTE: the pre-disinfection bath must be replaced after each operation.

### 9.2 Cleaning:

Automated cleaning process has been validated according to ISO 17664-1 and AAMI TIR 30.

The automated cleaning cycle to be used is as follows:

1. Pre-cleaning - temperature < 45°C - duration 2min minimum
2. Cleaning - temperature ≈ 55°C - duration 5min minimum
3. Neutralizing - approximate duration 2min
4. Rinsing with cold tap water - approximate duration 2min
5. Thermal Disinfection with deionized water - temperature ≈ 90°C - duration 5min minimum
6. Drying of the devices with a sterile gauze pad. Clean compressed air may be used to dry the lumens.
7. Perform a visual inspection on the devices and verify that they are clean, dry, and in proper working order prior to sterilization.

The automated cleaner need to be qualified according to the requirements defined in the ISO15883 standards.

Cleaning solution (detergent) has to be selected by the hospital to be compatible with the use in automated cleaner/disinfector according to the ISO 15883 standard and on surgical instruments made out of metallic and/or plastic parts.

Doses, temperatures and durations may be adapted according to solution manufacturer instructions.

Because of the use of some cleaning solutions, white traces may appear after drying.

A discoloration has no adverse effect on titanium alloy devices.

### 9.3 Inspection and functional test:

#### Instruments:

No maximal number of cycles is recommended but, between each use, all instruments need to be:

- Visually inspected for trace of organic residues corrosion (rust, pitting), damages such as scratches and notches, debris, discoloration, residue, flaking, wear, cracks and unreadable laser markings.
- Test for proper functioning of all joints, moving pieces and mechanisms and test the functionality of assemblies and connections. When needed, proceed to maintenance before sterilization.

Additional test for the Guiding wires: the straightness needs to be visually and functionally inspected (the guiding wire shall pass through the cannula of the screw).

Cannulated taps should be checked as bone material could obstruct the canula after drilling. A guiding wire should be inserted into the drill canula to eliminate any bone material and avoid obstruction.

Where appropriate, instruments are designed to be reassembled the identical way they were disassembled.

Any damaged instruments (for example with corrosion, excessive scratches, notches, residue, debris or unreadable laser markings) should be destroyed and replaced by a new instrument. In case of unexpected damage, please return the device to DISTIMP as a complaint.

Before sterilization check that the instruments are all there and are in good working order.

#### Implants and instruments:

For sterilization and storage, the devices can be settled in metallic containers and arranged in kits following carefully the DISTIMP recommendations to assure sterilization efficacy. Then, containers can be sterilized with the devices.

### 9.4 Sterilization:

All DISTIMP instruments and implants that are provided non-sterile and must be steam sterilized by hospitals and clinics according to one of the following validated methods.

| Method | Cycle      | Temperature   | Exposure time   | Drying time |
|--------|------------|---------------|-----------------|-------------|
| Steam  | Vacuum     | 134°C         | 18 min. minimum | -           |
| Steam  | Pre-vacuum | 270°F (132°C) | 4 min.          | 20 min.     |
| Steam  | Pre-vacuum | 134°C         | 3 min.          | 20 min.     |

Note for France: Use the first cycle parameters according to the French Circular DGS/ R13/2011/449 of 1 December 2011.

These methods are validated through ANSI/ AAMI ST79 to meet the requirements of a SAL of 10<sup>-6</sup>.

If after having followed this sterilization method there is still water in the sterilization containers or on/inside the device, the device must be dried, and sterilization repeated.

### 9.5 Maintenance:

The instruments joints, hinges, moving pieces and mechanisms should be regularly lubricated with a lubricant meant for medical grade surgical instruments, to reduce friction and wear.

Complying with the lubricant supplier's instructions is highly recommended. The lubricant used must be compatible with steam sterilization.

## 10. STORAGE AND DISPOSAL

- Store devices in a manner that provides protection from dust, light, insects, vermin and extreme temperature and humidity
- Avoid contact with chemical products or corrosive vapors.
- Do not store the instrument close to products that may have a corrosive action (chlorine).
- Unpacked devices can be stored in the dedicated container.
- It is IMPERATIVE to distinguish STERILE and NON-STERILE instruments.

The disposal of products shall be done according to the current procedures of the health care center to prevent any risk of cross-contamination.

## 11. WARRANTY

All warranty rights are lost if repairs or modifications are carried out by an unauthorized service center. The manufacturer does not take responsibility for any effects on safety, reliability or performance of the product if the product is not used in conformity with the instructions for use.

Technical alterations reserved.

## 12. COMPLAINTS/ INCIDENT AND INFORMATION

Any serious incident that occurs in relation to the device must be reported to the manufacturer as a complaint, and for Europe to the national competent authority of the patient/surgeons.

Any other grounds for dissatisfaction relating to the quality of the product should be notified to DISTIMP or its representative.

For all complaints, please give the name and reference along with the batch number of the component(s), your name and address and an exhaustive description of the event to help DISTIMP understand the cause of the complaint. When possible, please send the device to DISTIMP for the investigation.

All current information regarding VEOS products (instruction for use, surgical technique manual, check-list) is available on the website <https://www.distimp.com>.

DISTIMP SAS - 7 allée du Moulin Berger  
69130 ECULLY, FRANCE  
Tel : +33 (0)4 72 77 01 52  
[www.distimp.com](http://www.distimp.com)  
[quality@distimp.com](mailto:quality@distimp.com)

## SYMBOL KEY LEGEND FROM ISO 15223-1

Caution: symbol key is for reference only - some symbols listed may not apply. Please see main product label affixed to product.

MANUFACTURER

NON-STERILE

DO NOT REUSE

KEEP DRY

CATALOG NUMBER

KEEP AWAY FROM SUNLIGHT

CAUTION: CONSULT  
ACCOMPANYING DOCUMENTS

UNIQUE DEVICE IDENTIFIER

BATCH NUMBER

MR CONDITIONAL

CONSULT INSTRUCTIONS  
FOR USE

MEDICAL DEVICE

MANUFACTURING DATE

DO NOT USE IF PACKAGE IS  
DAMAGED AND CONSULT  
INSTRUCTIONS FOR USE

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