PROTHESE DISCALE / SANS CIMENT
 DISC PROSTHESIS/CEMENTLESS
 DISKUSPROTHESE / ZEMENTFREI
 PRÓTESIS DISCAL / SIN CEMENTO
 DISKUSPROTESE / CEMENTFRI
 VÄLILEVYPROTEESI / ILMAN SEMENTTIÄ
 PROTEZA DISKA / BEZ CEMENTA
 PROTESI DISCALE/SENZA CEMENTO
 SKIVEPROTESE / UTEN SEMENT
 DISCUSPROTHESE / ZONDER CEMENT
 PROTEZA KRAŻKA MIĘDZYKRĘGOWEGO / BEZCEMENTOWA
 PRÓTESE DISCAL / SEM CIMENTO
 PROTEZA DISCALĂ / FĂRĂ CIMENT
 DISKPROTES / CEMENTFRI
 PLATNIČKOVÁ PROTÉZA/BEZCEMENTOVÁ
 DISK PROTEZA / BREZCEMENTNA
 DISK PROTEZI / ÇİMENTOSUZ
 THIẾT BỊ ĐĨA ĐỆM/ KHÔNG CÓ XI MĂNG
 ΠΡΟΘΕΣΗ ΔΙΣΚΟΥ / ΧΩΡΙΣ ΤΣΙΜΕΝΤΟ
 人工椎间盘假体 / 非骨水泥型
 قرص بيففري اصطناعي / من دون مادة مُلصقة

LP.ESP

NOTICE D'INSTRUCTIONS
 INSTRUCTIONS FOR USE
 GEBRAUCHSANWEISUNG
 INSTRUCCIONES DE USO
 BRUGSANVISNING
 KÄYTTÖOHJEET
 UPUTE ZA UPORABU
 ISTRUZIONI PER L'USO
 BRUKERVEILEDNING
 GEBRUIKSAANWIJZING
 INSTRUKCJA STOSOWANIA
 FOLHETO DE INSTRUÇÕES
 MANUAL DE UTILIZARE
 BRUKSANVISNING
 NÁVOD NA POUŽITIE
 NAVODILA
 KULLANIM KILAVUZU
 HƯỚNG DẪN SỬ DỤNG
 ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
 说明书
 تعليمات الاستعمال

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WARNINGS AND INDICATIONS FOR USE
PLEASE READ THIS DOCUMENT CAREFULLY

Do not hesitate to consult our sales department if you need further information. This document is not exhaustive, nor is it a surgical operating manual presenting details about implantation. In all cases, please consult the included documents.

DESTINATION:

The devices covered in these instructions for use are intended to be implanted at the level of the lumbar spine between L3 and S1 during lumbar disc arthroplasty.

DESCRIPTION:

SPINE INNOVATIONS supplies the LP-ESP® lumbar disc prosthesis.

The LP-ESP® lumbar disc prosthesis is composed of two titanium alloy plates, separated by a deformable portion that is intended to allow the disc to move. The prosthesis has relief barbs and an HA (hydroxyapatite) coating on an undercoat layer of pure titanium (T40). The deformable portion is composed of a flexible polycarbonate urethane cushion and a central core made of silicon.

The LP-ESP® lumbar disc prosthesis is supplied sterile (gamma irradiation sterilisation) and is a single-use prosthesis.

Sales references	Type	Lordosis	Posterior height
255682	LP-ESP DISC PROSTHESIS / CEMENTLESS SIZE: 7° - H10mm	7	10
255683	LP-ESP DISC PROSTHESIS / CEMENTLESS SIZE: 7° - H12mm	7	12
255687	LP-ESP DISC PROSTHESIS / CEMENTLESS SIZE: 9° - H10mm	9	10
255688	LP-ESP DISC PROSTHESIS / CEMENTLESS SIZE: 9° - H12mm	9	12
255690	LP-ESP DISC PROSTHESIS / CEMENTLESS SIZE: 11° - H10mm	11	10
255691	LP-ESP DISC PROSTHESIS / CEMENTLESS SIZE: 11° - H12mm	11	12

IMPLANT COMPONENT MATERIALS:

The materials that the implants are made from are noted on each product label. The implant is radio-visible, except for the polycarbonate-urethane and the silicone. Qualitative composition of the implanted materials in mass percentage:

Sales references	Ti6Al4V	PCU	Silicone	T40	HAP
255682	63%<m<65%	25%<m<27%	7%<m<8%	1%	1%
255683					
255687					
255688					
255690					
255691					

PCU: Polycarbonate urethane

T40: Titanium

HAP: Hydroxyapatite

COMPATIBILITY BETWEEN DEVICES:

This is a monoblock device and does not require additional fixation elements.

PERFORMANCE / CLINICAL BENEFIT:

The LP-ESP® disc prosthesis must enable pain reduction, re-establish the lumbar curvature, and reproduce the functions of the disc.

The summary of safety and clinical performance (SSCP) of the LP-ESP® lumbar disc prosthesis is available in the European medical device database (Eudamed): <https://ec.europa.eu/tools/eudamed>.

The summary can be found on the Eudamed website using the UDI-DI of the LP-ESP® disc prosthesis: **376036923LP01BZ**.

INDICATIONS FOR USE:

The LP-ESP lumbar disc prosthesis is designed for the treatment of degenerative disc disease (DDD) in the lumbar column. Disc condition is evaluated based on the convergence of data from clinical examinations and imaging.

The LP-ESP® lumbar disc prosthesis is a medical device for primary treatment.

TARGET POPULATION:

Patients with a mature skeleton and less than 65 years of age failing to respond to at least 6 months of conservative treatment, presenting the aforementioned indications and presenting none of the contraindications listed below.

CONTRAINDICATIONS:

- Fractures
- Tumours
- Spinal stenosis, radiculopathy
- Significant segmental instability
- Spinal deformation, spondylolisthesis greater than 25%
- X-ray confirmation of severe lesions or degeneration of the facet joints
- Osteoporosis, osteochondrosis, and severe osteopenia
- Chronic or severe, local, spinal, or systemic infections
- Metabolic and systemic diseases
- Sensitivity to the implant's materials
- Dependence on drugs, drug addiction or alcoholism
- Pregnancy
- Obesity
- Lack of cooperation from the patient
- Pathologies and surgical situations that preclude any benefit of spinal surgery, such as arteritis of the lower limbs, genitourinary disorders of neurological origin, damage to more than 2 lumbar discs.

PRECAUTIONS:

A - GENERAL PRECAUTIONS

- The lifespan of implants is affected by numerous biological, biomechanical, and other factors. As a result, carefully following the indications, contraindications, precautions, and warnings for this product plays an essential role in its use.

The outcome of intervertebral disc prosthesis depends on the patient's medical history. Psychological preparation is indispensable.

Patients must be informed of the limitations of the device, including among others, the impact of overburdening caused by weight or excessive activities. They must be advised as to how to modify their activities accordingly. Joint replacement can never reproduce the functions formerly performed by a normal, healthy joint. The patient should not hold on to unrealistic hopes about its functionality and should consult the surgeon in the event of abnormal problems experienced at the device loca-

tion.

- Before clinical use, the surgeon and the surgical team must be trained in using the device and its ancillary equipment. They must have assimilated the aspects of the surgery, as well as the limitations of the device.
- Care should be taken to ensure that the components and all polished or coated surfaces are protected against abrasion, scratches, or any other damaging effects caused by abrasive or metallic objects.
- The contact zones of the implant must be clean and free of debris (bone, cement) before implantation. The surgeon should make certain to obtain even contact of the implant on the bone.
- It is advisable to make sure that the presence of other devices does not compromise the device's functional integrity.
- Implants provided by SPINE INNOVATIONS must not be mixed with implants of another origin, or with different systems.

B - CEMENTLESS IMPLANTS

During handling, great care should be given to the coating. Any contact with the coating before implantation must be avoided. Biologically fixed implants (identification on the product label: T40 + HA) must not be cemented.

IMPLANTATION TECHNIQUE:

The specific ancillary equipment, delivered non-sterile, is necessary for implanting implants designed by SPINE INNOVATIONS. The surgical techniques are available from SPINE INNOVATIONS or its distributors.

If hesitating between two sizes, use the smaller prosthesis if the test prosthesis is sufficiently maintained between the vertebrae (see surgical technique).

CONDITIONS OF USE:

The prostheses must be implanted in an operating facility, under aseptic conditions, and in compliance with hygienic practice, by a surgeon specialised in the treatment of spinal pathologies who regularly practices this type of implantation. Implantation must be carried out using the appropriate ancillary equipment, with undamaged instruments, following the recommended surgical technique and the instructions for use, to treat a patient who presents with the indications defined above.

WARNINGS:

Patients receiving an implant must be informed that the implant's lifespan may depend on their weight and activity level. Damaged or defective implants must not be used.

The design of this device does not allow for its re-use; re-use may alter its mechanical characteristics, which may affect its performance. Re-use of a single-use device may also result in cross-contamination, thus causing patient infection.

Using an implant that is of the wrong size or sacral tilt runs the risk of decreasing its stress resistance. SPINE INNOVATIONS is not liable for use of another manufacturer's components, or use of its own components for purposes other than those intended.

The implants must not be modified or subjected to any treatment.

ADVERSE EFFECTS:

It is the surgeon's responsibility to provide the patient with all information before the operation, and especially to inform the patient of the following potential adverse effects:

- Sexual disorders
- Pain
- Degeneration of adjacent discs
- Implant migration
- Accumulation of liquid or formation of a hernia in the operated area
- Infection of the wound and/or systemic infection
- General reoperation
- Blood circulation disorders
- Surgical error

- Blood clots and blood flow restriction liable to cause a pulmonary embolism
- Trauma caused by the surgical intervention, such as nerve damage or damage to the spinal cord, excessive bleeding and/or fracture of the vertebral body
- Reoperation for implant revision
- Reoperation for additional fixation
- Reoperation for implant removal
- Complication of the bone and/or tissue structure (spondylolisthesis, degeneration of the facets, herniation of an adjacent disc, spinal stenosis)
- Sagging of device
- Neurological disorders
- Device breakage
- Nerve damage or damage to the spinal cord liable to cause a loss or deficiency of sensory function
- Cardiovascular problems liable to cause a heart attack or stroke
- Ossification
- Periprosthetic fracture
- Allergy
- Osteolysis

IMPORTANT INFORMATION TO GIVE THE PATIENT:

The patient must be provided with:

- Information to enable device identification, including the name, batch number, UDI, device model, as well as the name, address, and website of the manufacturer.
- Warnings, precautions, and measures to be taken by the patient or by healthcare professionals regarding mutual interference with external environmental sources or conditions or reasonably predictable medical examinations.
- Information on the expected lifespan of the device and necessary monitoring is provided on the patient instructions, accessible using the QR code on the implant card given to the patient.
- Further information intended to guarantee safe use of the device by the patient is provided on the patient instructions, accessible using the QR code on the implant card given to the patient.

This information is indicated on the implant card and on the patient instructions, accessible using the QR code on the implant card.

Access to patient instructions is provided via the QR code below or at the following address: <https://spine-innovations.com/patient-information>

MRI COMPATIBILITY:

Non-clinical tests have shown that the ESP disc range is "MR conditional", according to the definitions of norm ASTM F2503-20. A patient with a device from this range can therefore have a totally safe scan using an MRI system that meets the following conditions:

- Patient implanted with a single ESP disc
- Patient with no thermoregulation disorder (i.e. no alteration of systemic thermoregulation or no reduction of local thermoregulation) and patient in controlled conditions (one specific doctor or trained person able to respond immediately to the physiological stress induced by the heat).
- MRI system with horizontal tunnel and static magnetic field of 1.5 Tesla or 3 Tesla.
- Spatial gradient of the magnetic field less than or equal to 19T / m.
- Product B0 * IdB 0 / drl less than or equal to 48T² / m.
- Use of the whole-body emission/reception RF coil only.
- First level controlled operating mode, i.e. SAR (specific absorption rate) averaged for the whole body (WB-SAR) limited to 4 W / kg.
- During non-clinical tests, after 15 minutes of continuous acquisition, the ESP disc produced a maxi-

mal temperature increase of 5.0 ± 1.0 °C at 1.5 T for measured WB-SAR of 3.50 ± 0.81 W / kg and a maximal temperature increase of 3.5 ± 1.0 °C at 3 T for measured WB-SAR of 3.94 ± 0.88 W / kg.

- For WB-SAR of 4W / kg, the ESP disk can be expected to produce a maximal temperature increase of 5.7 ± 1.8 °C at 1.5T and 3.6 ± 1.3 °C at 3T.

The quality of MRI images may be compromised if the area of interest for the imaging is in the same zone as the implant.

It may be necessary to adjust the acquisition parameters to compensate.

IMPORTANT INFORMATION FOR THE SURGEON:

Any serious accident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user or the patient has residence.

PACKAGING AND STERILISATION:

The information on the product label ensures manufacturing traceability and indicates which sterilisation method has been used.

IMPLANTS DELIVERED STERILE:

Verify that the packaging elements (peelable pouches, blister packs and films) are perfectly sealed and that all elements are undamaged before using the implants. Do not use a product if its packaging is damaged or if the tamper-proof seal has been broken. The product must not be resterilised by any method whatsoever at the risk of altering its mechanical characteristics.

The sterilisation indicator sticker on the outer packaging should be red for radiation sterilisation. In all cases, an indicator showing the pre-sterilisation colour, yellow, may indicate that the product is not sterile, and, in this case, the product must be returned. Desterilised instruments cannot be returned.

The expiry date is noted on the product label.

DEVICE DISPOSAL:

In the event of explantation, the implant must be sent to a specialised waste handling service so that it can be disposed of properly with respect to the environment and strict sanitation laws. A defective device that has been removed must be decontaminated and returned to the manufacturer.

DATE OF FIRST CE MARKING:

LP-ESP disc prosthesis: 12-2023

The reference text is the French text.

Symbols meaning	
	Distributor
	Manufacturer
	Consult instructions for use
	Caution
	Sterilised using irradiation
	Protect against humidity
	Keep away from sunlight
	Do not use if packaging is damaged
	Do not reuse
	Do not resterilise

	MR Conditional
Rx Only	Prescription only
	Date of manufacture
	Use by date
	Batch code
	Catalogue reference
	Double sterile barrier system
	Unique device identifier
	Medical device
	Patient name or identification number
	Implantation date
	Healthcare facility name and address
	Patient information website

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